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EN ATTENTION!!!

This is a single use product. Reutilization is not recommended. Do not use if the packaging is open or damaged. Non-sterile product. Must be sterilized before use. Single use recommended by the manufacturer.

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1.1 Modelo comercial: "Tampa Pilar Hexagonal"

Commercial model: "Hexagonal Abutment Healing Cap"

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EN Product description

The Hexagonal Abutment Healing Cap is a straight abutment made from titanium F67. The Hexagonal Abutment Healing Cap is a single-use piece. It is compatible with different Hexagonal Abutment models, and is used as a protection during the healing period, allowing the gum tissue to acquire the best shape for the desired aesthetic profile.

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This product is individually packaged. Each package contains:

a) Titanium F67 Healing cap and three medical chart labels

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EN Specifications and technical features

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EN Raw Material: The raw materials used for the "Hexagonal Abutment Healing Cap" is commercially pure titanium ASTM F-67

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Commercial model "Mini Tapered Healing Cap"

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EN Description

The Mini Tapered Healing Cap is a straight abutment made from titanium F67. The Mini Tapered Healing Cap is a single-use piece. Its base has the shape of a 1.17 hexagon for connection with the insertion tool and a M1.4 screw to fix the healing cap to the corresponding component. It is compatible with different Mini Tapered Abutment models, and is used as a protection during the healing period, allowing the gum tissue to acquire the best shape for the desired aesthetic profile.

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EN Raw Material: The raw materials used for the "Mini Tapered Healing Cap" is commercially pure titanium ASTM F-67

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1.3 Modelo comercial "Tampa Cone Estético"

EN Commercial model "Tapered Aesthetic Healing Cap"

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EN Description

The Tapered Aesthetic Healing Cap is a straight abutment made from titanium F67. The Tapered Aesthetic Healing Cap is a single-use piece. Its base has the shape of a 1.17 hexagon for connection with the insertion tool and a M1.4 screw to fix the healing cap to the corresponding component. It is compatible with different Tapered Aesthetic Abutment models, and is used as a protection during the healing period, allowing the gum tissue to acquire the best shape for the desired aesthetic profile.

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
EN Accompanying parts and accessories

This product is individually packaged. Each package contains a Titanium F67 Healing cap and three medical chart labels

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EN Raw Material: The raw materials used for the "Tapered Aesthetic Healing Cap" is commercially pure titanium ASTM F-67

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EN Commercial Model "Cone Morse Trunion Cover"

A) Description

The Cone Morse Trunion Healing Cover is comprised of a straight pillar, made of Titanium F67. The Cone Morse Trunion Healing Cover is for single use. In its upper base, there is a 1.17 hexagon for connection to the placing key. In its lower base, there is a M1.4 thread for screwing to the corresponding component. It is compatible to the different models of Cone Morse Trunion Prosthetic Pillar and protects it during the healing period, in order to allow the gingival tissue to obtain a shape adequate to the desired esthetic profile.

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EN B) Parts and accessories included

The product is packed individually

Each package is composed of:

- one Healing Cover made of Titanium F67

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EN C) Specifications and technical characteristic

Dimensions	Base Ø (mm)	Belt height (mm)
Cone Morse Trunnion Cover	2.8; 3.5; 4.5	4.0 and 6.0

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EN Raw Material:

Commercially pure Titanium ASTM F-67 is used as raw material for the "Cone Morse Trunnion Cover".

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EN Physical principle and product technology fundamentals, applied to their operation and action indication:
This product works as protection for the areas to be used as link between the prosthetic component and the implant itself.

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EN Accompanying parts and accessories:

This product is individually packaged. Each package contains 1 Implacil De Bortoli Healing Cap and 3 Medical chart labels.

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EN Optional, consumables and support materials

Optionals: Not applicable. Consumables: Not applicable. Items manufactured to be used exclusively with the company's products: All pieces described in these Usage Instructions are to be used exclusively on Implacil Healing Caps.

Attention: The use of any piece, accessory or material not specified or described in these Usage Instructions is of the entire responsibility of the user.

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EN Directions for product usage

1. The Implacil De Bortoli Healing Cap is packaged inside a Tyvek sealed blister.

2. Check the packaging to see if it is in perfect conditions. Do not use the product if packaging is damaged.

3. Remove the Healing Cap from the blister.

4. The Healing Cap is pre-sterilized in a classified room. To sanitize it again before its use, follow the steps below: rinse the item with cold water; wash it with a neutral to basic detergent; rinse it again with deionized water or water purified by reverse osmosis.

This procedure aims to completely eliminate microorganisms such as viruses, bacteria, microbes, fungi (vegetative or sporulating). After the drying process, the instruments should be sterilized by autoclaving, using temperature and exposition time according to the following method: Sterilization temperature: 132°C, Sterilization cycle: 20 minutes, and Drying cycle: 20 – 60 minutes. The surgical instruments should be placed inside the container (with holes on the lid and the bottom), or in a tray or box; as long as they are placed as to be prepared for the sterilization process. Do NOT open the autoclave to avoid rapid condensation; let the vapor out completely so that the Drying cycle is completed by itself. Clean the autoclave thoroughly and periodically, removing all dirt and rust residue formed (iron oxide).

5. The product is non-sterile, and should be sterilized by the professional according to these steps:

Sterilization temperature: 132°C, Sterilization cycle: 20 minutes and Drying cycle: 20 – 60 minutes. The items should be placed inside the container (with holes on the lid and the bottom), or in a tray or box; as long as they are placed as to be prepared for the sterilization process. Do NOT open the autoclave to avoid rapid condensation; let the vapor out completely so that the Drying cycle is completed by itself. Clean the autoclave thoroughly and periodically.

6. Screw the Healing Cap to the component to be protected.

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EN Special conditions for product transportation, storage and packaging

Storage: Keep away from direct contact with sun and rain, and inside its original packaging. Transportation: protect the packaging from rain. Packaging: Keep away from direct contact with sun and rain.

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EN Warnings and/or precautions to be adopted

Storage: Keep the product inside its original packaging and away from direct contact with sun and rain. This product is not heat-sensitive and exhibits no variation in environments with high or low levels of humidity, which is why it can be sterilized by autoclave. Even so, it is recommended that you keep the product between 15°C and 35°C, with humidity between 45% and 80%, to ensure the packaging's stability. **Transportation:** The Healing Caps must be transported in clean and dry conditions at room temperature. Transportation must be adequate in order to avoid falls or damages in the original packaging. **Packaging:** Implacil De Bortoli Healing Caps can only be used by specialized professionals in the area. It is necessary to elaborate a prosthetic execution plan, observing the right diameter, angle and length for the adaptation of the prosthetic component to the implant in order to avoid future complications.

Implacil De Bortoli Healing Caps were designed for single use, and their reprocessing is prohibited. Do not use it in any way if the packaging is violated. The product is not sterile, so a sterilization procedure must be performed by the professional before any intervention on the patient.

The professional must observe the pressure exerted in the moment of application so as to not damage the item, and inform the patient about the necessity of periodic monitoring and avoiding physical and mechanical efforts in order not to submit the product to undue efforts, as well as the necessity of proper hygiene.

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EN 5 - Performance on essential medical product safety and efficacy requirements and eventual undesirable side effects

5.1 - Indication, purpose or use for the product

Implacil De Bortoli Healing Caps are indicated for installation over the abutments, so as to offer a transition from the installation of the prosthetic components to the preparation of the implant. It protects the areas to be adapted, acting as a mechanical barrier to gum tissue growth and allowing access to the abutment's platform.

5.2 - Undesirable secondary or side effect and contraindications

The product application final result may be unsuccessful if an inadequate technique is used or if the patient is not subjected to the proper examinations and prosthetic planning. In these cases the product might be lost or fractured. Incorrect application of the product may cause side effects in the place where it was applied, such as pain, swelling, sensitivity, tissue reaction and infection.

Implacil De Bortoli Healing Caps are not indicated in cases of patients with chronic periodontal inflammation, inadequate oral hygiene or parafunctional habits, active intra-oral infection or in cases of immediate load where the primary stability of the implant is not adequate.

5.3 - Product safety and efficacy:

Implacil De Bortoli Healing are completely safe as long as the safety rules are observed in all procedures performed. The product was developed in non-toxic material for odontologic use, manufactured in Titanium F67. The product's physical and dimensional characteristics will not change as long as it's stored and kept according to its usage indications. The product is not sterile, so a sterilization procedure must be performed by the professional before any intervention on the patient.

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EN 6 – Installation or connection to other products

Implacil De Bortoli Mini Tapered Healing Cap and Tapered Aesthetic Healing Cap must be connected to an insertion tool and screwed to the component to be protected. Implacil De Bortoli Hexagonal Abutment Healing Cap must be screwed to the component to be protected.

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EN 7 – Concerns in case of product disposal

Para evitar contaminación ambiental o improprio uso de Implacil Healing Caps, the ones not used must be discarded in any specific location for potentially contaminated materials.

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EN 8 – Term of Warranty

IMPLACIL DE BORTOLI MATERIAL ODONTOLOGICO LTDA. assures the owner of this product warranty against any material or manufacture defects, if presented in a period of 365 days, from the date of acquisition by the first acquirer.

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